(en)

air-Q SP 3/3G Placement Procedure

The procedure below is intended as a guide. Many techniques can be successfully used to properly place the air-Q SP 3/3G into the pharynx.

- Lubricate the external surface including the mask cavity ridges.
- Open the patient's mouth and elevate the tongue. Elevating the tongue lifts the epiglottis off the posterior pharyngeal wall and allows the air-Q SP 3/3G easy passage into the pharynx. A mandibular lift is especially recommended. A tongue blade placed at the base of the tongue also works well for this purpose.
- Place the front portion of the air-Q SP 3/3G mask between the base of the tongue and the soft palate at a slight forward angle, if possible.
- Pass the air-Q SP 3/3G into position within the pharynx by gently applying inward and downward pressure, using the curvature of the air-Q SP 3/3G mask and airway tube as a guide. Simply rotate the air-Q SP 3/3G forward and inward. Minimal manipulation may be necessary to turn the corner into the upper pharynx. Continue to advance until fixed resistance to forward movement is felt. Correct placement is determined by this resistance to further advancement. Some users place the back of the left index finger behind the mask, flexing the finger forward to help guide the mask around the corner into the pharynx. Once the mask has negotiated the turn, the left hand is then used to do a mandibular lift while exerting downward and inward pressure on the air-Q SP 3/3G with the right hand during final advancement into the pharynx. This technique seems to be easy to learn and is particularly successful.
- Tape the air-Q SP 3/3G in place.
- Check the air-Q SP 3/3G connector to ensure it is fully engaged within the airway tube, and attach the connector to the appropriate breathing device. Check for adequate ventilation.
- Place a bite block between the patient's teeth. Keep the bite block in place until the air-Q SP 3/3G is removed.

air-Q SP 3/3G Intubation Procedure

The air-Q SP 3/3G by Cookgas® LLC is intended not only to be an outstanding airway for general use, but also to be a simple and reliable tool for intubation of the trachea with OETTs. Due to its patented design, standard OETTs (sizes 9.0 mm - 3.0 mm) can be easily passed through the air-Q SP 3/3G and into the trachea. Further, the air-Q SP 3/3G can be easily removed following intubation with the aid of the patented air-Q Removal Stylet, also by Cookgas® LLC. The following procedure for intubation is only intended as a guide. Many techniques can be successfully used for tracheal intubation using the air-Q SP 3/3G.

- Prior to intubation, the laryngeal musculature and vocal cords must be relaxed, either by an aerosolized local anesthetic or with the aid of a muscle relaxant.
- Pre-oxygenate.
- Prepare the appropriately-sized OETT by completely deflating the OETT cuff and lubricate well. It is important to deflate the OETT cuff completely to allow the OETT to slide easily within the air-Q SP 3/3G.

air-Q SP 3/3G Removal Procedure

Removing the air-Q SP 3/3G following OETT intubation is easily accomplished with the aid of the air-Q Removal Stylet by Cookgas® LLC. The air-Q Removal Stylet consists of an adapter connected to a rod. The adapter is tapered from bottom to top, and has horizontal ridges and vertical grooves. The taper allows the stylet to accommodate multiple OETT sizes. The ridges engage the OETT in a firm, secure grip, giving the user control of the OETT during the air-Q SP 3/3G removal process. The grooves allow spontaneously breathing patients unimpeded air passage within the OETT during removal of the air-Q SP 3/3G. By immobilizing and exerting an inward stabilizing force on the OETT, the air-Q SP 3/3G Removal Stylet allows for the swift removal of the air-Q SP 3/3G without dislodging the previously-placed OETT from the patient.

- Remove the OETT connector from the OETT.
- Squeeze the proximal portion of the OETT between the index finger and the thumb, leaving enough room for the adapter portion of the stylet to enter the proximal opening of the OETT. Alternatively, squeeze the proximal end of the air-Q SP 3/3G airway tube, trapping the OETT inside.
- Insert the tapered end of the air-Q Removal Stylet into the proximal OETT (the long axis should be in the 12 o'clock - 6 o'clock position) until the adapter fits snugly within the OETT.
- For larger sizes, (2.0 - 4.5), with firm inward pressure, rotate the stylet adapter in a clockwise direction (into the 3 o'clock - 9 o'clock position) until the adapter firmly engages the OETT. For smaller sizes (0 - 1.5) simply push the stylet firmly into the OETT. Please practice this a few times prior to attempting on a patient.
- Deflate and lubricate the pilot balloon on the OETT prior to withdrawing the air-Q SP 3/3G. Reinflate the OETT following air-Q SP 3/3G removal.
- While exerting an inward stabilizing force on the stylet, slowly withdraw the air-Q SP 3/3G outward over the stylet rod.
- For larger sizes (2.0 - 4.5) pass the stylet through and through. For smaller sizes (0 - 1.5) remove the stylet from the proximal end of the air-Q SP 3/3G while stabilizing the OETT at the mouth. Discard single-use air-Q SP 3/3Gs following use.
- Reposition the OETT to the proper depth within the patient, if needed, and then tape into place.
- Replace the OETT connector within the OETT. Inflate the OETT, if needed, and attach to an appropriate breathing device. Check for adequate ventilation.

Cautions/Warnings

- Inspect all air-Q SP 3/3G devices prior to use. Discard defective devices.
- Do not use sharp instruments on or near the air-Q SP 3/3G.
- Confirm that the air-Q SP 3/3G size matches the connector size prior to use.
- The air-Q SP 3/3G connector is removable and airway disconnect is possible. Take standard precautions to minimize the chance of disconnect.

Warranties

Cookgas® LLC agrees to warrant the disposable air-Q SP 3/3G for a period of 30 days following the invoice date. Warranty covers materials and manufacturing defects provided that the airway is used according to the procedures outlined in the Instructions For Use (IFU) manual. Warranty is valid only following purchase from authorized distributors.

Cookgas® LLC disclaims all other warranties whether expressed or implied.

Patents U.S.	Patents CAN	Patents UK
5,937,860	US 6,705,321 B2	2,231,331
US 6,422,239 B1	US 7,357,845 B2	2,371,435
US 6,892,731 B2	US 7,780,900 B2	2,410,043
US 7,331,347 B2	US 7,784,464 B2	
US 7,900,632 B2	US 7,934,502 B2	Patents JAPAN
US 9,320,864 B2	US 8,622,060 B2	5296058
	US 8,631,796 B2	
	US 10,729,866 B2	

Other USA & Foreign Patents Pending

ATTENTION
SEE
DIRECTIONS
FOR USE

 Latex Free

 Keep Dry

 Keep
Away From
Sunlight

Rx
Only
By
Prescription
Only

 EXP
Expiration
Date

 Single Use
DO NOT
REUSE

 MR Safe

 Do Not
Resterilize

Distributed exclusively by:

For Ordering Information, Contact:

2710 Northridge Dr. NW, Ste A
Grand Rapids, MI 49544 USA
PH. 800-433-2797
www.Sun-Med.com

M T Promedt Consulting
Altenhofstrasse 80
D-66386 St. Ingbert • Germany
Tel. +49(0) 6894.581020 • Fax +49(0).6894.581021
email: info@mt-procons.com
www.mt-procons.com

The air-Q SP 3/3G represents a revolution in airway device design. It functions similarly to the original air-Q in that it retains the soft perimeter mask cuff. This gives the mask the ability to change its size and shape depending on each patient's pharyngeal anatomy, improving fit. Like the original, it can also be used as a conduit for intubation in difficult airway situations. It differs in that it does not contain an inflation valve and so does not inflate normally. Instead, the air-Q SP 3/3G incorporates a large aperture between the internal volumes of the breathing tube and the peripheral cuff. This aperture allows fluid communication between the tube and cuff to "self-pressurize", causing inflation of the cuff during positive pressure ventilation.

This product is to be used by trained personnel only.

Available in Single Use Only


Recommendations:

Size	IBW	Max. OETT	Mouth Opening†	← →	Max OG Tube‡
5	≥80-kg	9.0 mm	25 mm	20 cm	18
4	60-80 kg	8.0 mm	23 mm	18 cm	16
3	30-60 kg	7.0 mm	20 mm	16 cm	14
2	17-30 kg	5.5 mm	17 mm	14 cm	10
1.5	7-17 kg	5.0 mm	14 mm	11 cm	8
1.0	4-7 kg	4.5 mm	11 mm	9 cm	8
0.5	2-4 kg	4.0 mm	8 mm	7 cm	6
0	< 2 kg	3.0 mm	5 cm	6 cm	5

- Minimum mouth opening for insertion.
- Distance from the external edge of the airway tube to the internal ventilatory opening.
- Maximum OG size (air-Q SP 3G only).

Recommended insertion Technique

Anbefalet indføringsteknik
Aanbevolen inbrengtechniek
Technique d'insertion recommandée
Empfohlenes Einführungsverfahren
Συμμετώμενη τεχνική εισαγωγής
Tecnica di inserimento consigliata
Anbefalt innføringsteknikk
Zalecana technika wsuwania
Técnica de Inserção Recomendada
Рекомендованная техника введения
Técnica de inserción recomendada
Rekommenderad införselsteknik
Önerilen Sokma Tekniği
推荐插入方法


Recommended Depth of Insertion Range

Anbefalet område for indføringsdybde
Aanbevolen inbrengdiepte bereik
Profondeur recommandée de la portée d'insertion
Empfohlener Einführtiefenbereich
Συμμετώμενο βάθος εύρους εισαγωγής
Gamma di profondità di inserimento raccomandata
Anbefalt område for innføringsdybde
Zalecany zakres głębokość wsuwania
Intervalo de Profundidade de Inserção Recomendada
Рекомендованный диапазон глубины введения
Profundidad recomendada del rango de inserción
Rekommenderat djup för införselsspänn
Önerilen Sokma Aralığı Derinliği
推荐插入深度的范围

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