Instructions for Use: Endotracheal Tube

Not Made with Natural Rubber Latex

EN Endotracheal Tube

Federal (USA) Law restricts this device to sale by or on the order of a physician or any other practitioner licensed by the law of the State in which he/she practices to use or order the use of the device.

Intended Use: Endotracheal Tubes (cuffed or un-cuffed) are indicated airway management.

Description:

- Endotracheal Tube (cuffed or un-cuffed)
- Smooth Murphy Eye to provide an alternate ventilation port if tube becomes occluded.
- Softer rounded Murphy Eye is less invasive.
- Smooth tip to minimize trauma during intubation.
- Radiopaque strip embedded to assess exact location of tube.
- High volume, low pressure cuff provides an effective low pressure seal.
- Standard 15mm connector supplied with each tube, which insures compatibility with circuit connectors.
- Bold, easy-to-read depth marks for placement in low light conditions.
- Valve with sensitive pilot balloon indicates cuff's status.

Material: Non-irritating, Non-sensitizing, Non-cytotoxic

Contraindications: As with any PVC Endotracheal Tubes, use of Endotracheal Tubes in procedures which will involve the use of a laser beam or electrosurgical active electrode in the immediate area of the device is contraindicated.

Warnings/Cautions:

- Verify proper assembly/attachment.
- Test cuff by inflating it prior to use.
- Deflate cuff prior to insertion.
- Tube placement and intra-cuff pressure should be continuously monitored.
- A second method of monitoring tube placement is recommended.
- Safe cuff pressure should be less than or equal to 25 cm H₂O.
- Reassess tube placement whenever patient movement, clinical interventions, or clinical condition suggest the possibility of tracheal tube dislodgement.
- Beware of breathing circuit disengagement when transferring patients by stretchers.
- Use proper tube fixation methods.
- Do not use lubricants with a petrochemical base.
- Single use only.
- Do not resterilize.) Resterilization could change the physical characteristics of the material and should not be attempted.
- Do not use if package has been previously opened or damaged.
- Exposure to elevated temperatures and ultraviolet light should be avoided during storage.
- Store at a temperature that does not exceed 49°C.
- Sterile if package is not opened or damaged.

Suggested Directions for Use:

- Remove the sterile device from the package and retain package for patient record.
- Test the cuff inflation system prior to intubation (in the case of cuffed Endotracheal Tube).
- Deflate cuff prior to intubation (in the case of cuffed Endotracheal Tube).
- If using a stylette, verify that the end of the stylette does not protrude from the end of the tube or pass through the Murphy Eye.
- Inflate the cuff slowly while using the minimum amount of air required to provide an effective seal (in the case of cuffed Endotracheal Tube).
- Monitor the cuff pressure levels during intubation in order to insure the effective seal (in the case of cuffed Endotracheal Tube).
- Tracheal tube placement should be confirmed by viewing the position of the tube tip with chest radiograph.
- Remove the tube only after complete deflation of the cuff (in the case of cuffed Endotracheal Tube).
- Dispose device and device packaging in accordance to hospital practices and applicable state, federal and other regulations.

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SunMed, LLC
Manufactured for Ventlab, LLC,
a Subsidiary of SunMed Holdings, LLC
Mdi Europa GmbH | Langenhagener Str. 71