




AIR-Q®SP3G SIZING

ITEM	SIZE	IBW (LBS)	IBW (KG)	MAXIMUM OETT	MOUTH OPENING	LENGTH OF TUBE	MAXIMUM OG TUBE	AIR Q3 CONNECTOR
60005	0	< 4 lbs	< 2 kg	3.0 mm	5 mm	6 cm	5 Fr	
60055	0.5	4 - 9 lbs	2 - 4 kg	4.0 mm	8 mm	7 cm	6 Fr	
60105	1	9 - 15 lbs	4 - 7 kg	4.5 mm	11 mm	9 cm	8 Fr	
60155	1.5	15 - 37 lbs	7 - 17 kg	5.0 mm	14 mm	11 cm	8 Fr	
60205	2	37 - 66 lbs	17 - 30 kg	5.5 mm	17 mm	14 cm	10 Fr	
60305	3	66 - 132 lbs	30 - 60 kg	7.0 mm	20 mm	16 cm	14 Fr	
60405	4	132 - 176 lbs	60 - 80 kg	8.0 mm	23 mm	18 cm	16 Fr	
60505	5	> 176 lbs	> 80 kg	9.0 mm	25 mm	20 cm	18 Fr	