


HYPERINFLATION SYSTEM

LATEX FREE | SINGLE PATIENT USE | DISPOSABLE

The Hyperinflation System is available in five different pleated inflation bag sizes. Available with Ventlab's signature integrated pressure manometer and components to create a custom system.



HYPERINFLATION SYSTEM

- System equipped with a stay-put dial to enable clinician to achieve targeted static pressure
- Five pleated breathing bag sizes available: 1/4, 1/2, 1, 2, 3 liter
- Seven inflatable face mask sizes available: small infant, infant, small child, child, adult (regular, large, wide)
- Non-inflatable face masks can be specified

ITEM	DESCRIPTION	PK
HS1000	Hyperinflation System	20
HS4000	Hyperinflation System + Integrated Manometer	20

HYPERINFLATION SYSTEM

CUSTOMIZE YOUR RESUSCITATION / VENTILATION BAG


Choose from the following components to creat a custom hyperinflation system to meet your needs



Hyperinflation
4000-Series

Shown with integrated
pressure manometer
and pop-off valve

SAFETY COMPONENTS



INTEGRATED PRESSURE MANOMETER

Ventlab's innovative color-coded dial clearly denotes pressure delivery:

- Green - Target pressure level <20 cm H₂O
- Yellow - Potentially dangerous pressure level
- Red - Dangerous pressure level



POP-OFF VALVES

40 cm H₂O & 25 cm H₂O
Allows clinician to use an appropriate pop-off level

CONVENIENCE COMPONENTS



FACE MASKS

Inflatable and non-inflatable
face masks available