

PREMIUM ERGO MASK

LATEX FREE | SINGLE PATIENT USE | DISPOSABLE



No air cushion design means no need to inflate; no risk of puncture or leaking



Soft, flexible seal designed to fit the face more precisely, more comfortably

Curved mask shape mimics facial landscape for a better patient fit



Ergonomic design allows for comfortable handling and control for clinician

Finger grips for secure mask positioning

PREMIUM ERGO MASK

- Ergonomic shape combined with finger grips on both front and back of mask provide a secure hold when positioning
- Anti-slip finger grips feature help maintain secure hold even when hands are wet
- Unique anatomical shape and seal enables a firm, yet comfortable patient fit
- Bright color-coding for accurate size selection
- Specifically designed without a cushion to minimize deadspace
- Can be used at higher elevations with non-air cushion design
- Airtight fit available in 6 sizes



ITEM	DESCRIPTION	COLOR	SIZE	PK
VR0900	ERGO Mask	White	Small Neonate	20, 40
VR1900	ERGO Mask	Gray	Small Infant	20, 40
VR2900	ERGO Mask	Orange	Infant	20, 40
VR3900	ERGO Mask	Yellow	Small Child	20, 40
VR4900	ERGO Mask	Blue	Child	20, 40
VR5900	ERGO Mask	Red	Adult	20, 50