


# ETCO<sub>2</sub> CANNULA + SCOOP


LATEX FREE | SINGLE PATIENT USE | DISPOSABLE | CE MARK

End Tidal CO<sub>2</sub> Cannula simultaneously delivers oxygen while capturing CO<sub>2</sub>. A series of strategically positioned holes deliver O<sub>2</sub>, and CO<sub>2</sub> is effectively captured through the nares and scoop. This design helps lessen attenuation of the end tidal sample by oxygen dilution, particularly for mouth-breathing patients.



## ETCO<sub>2</sub> CANNULA + SCOOP

- Innovative design suitable for mouth-breathing patients
- Dual nasal cannula simultaneously provides oxygen to both nostrils while obtaining CO<sub>2</sub> sampling during spontaneous breathing
- Designed to prevent mixing of fresh oxygen with CO<sub>2</sub>
- 3-Channel tubing helps prevent kinking
- Includes male/female adapter with luer lock connector
- Available in three tubing lengths



ITEM	DESCRIPTION	SIZE	LENGTH	PK
4707F	EtCO <sub>2</sub> Cannula + Scoop & Fits-All Connector	Adult	7'	25
4710F	EtCO <sub>2</sub> Cannula + Scoop & Fits-All Connector	Adult	10'	25
4714F	EtCO <sub>2</sub> Cannula + Scoop & Fits-All Connector	Adult	14'	25

ITEM	DESCRIPTION	SIZE	LENGTH	PK
4707	EtCO <sub>2</sub> Cannula + Scoop	Adult	7'	25
4710	EtCO <sub>2</sub> Cannula + Scoop	Adult	10'	25
4714	EtCO <sub>2</sub> Cannula + Scoop	Adult	14'	25