


DUAL LUMEN ENDOBRONCHIAL TUBE

LATEX FREE | STERILE | SINGLE USE | DISPOSABLE

The Robertshaw dual lumen endobronchial tube is available in left and right designs. Helps provide an airtight seal of the trachea and one bronchus. Packaged with dual connectors, suction ports, suction catheter, comes sterile, and fully assembled.



DUAL LUMEN ENDOBRONCHIAL TUBE

- Color coded clear (tracheal) and blue (bronchial) for easy differentiation of appropriate lumen, pilot balloon and cuff
- Soft PVC
- High volume, low pressure cuff
- Radiopaque strip
- Graduated depth marks
- Suction catheter included
- Bronchial stylette included
- 15 mm Y-connector included with each tube

DUAL LUMEN (LEFT)

ITEM	DESCRIPTION	SIZE	PK
1-7370-35	Endobronchial Robertshaw Left	35 Fr	1
1-7370-37	Endobronchial Robertshaw Left	37 Fr	1
1-7370-39	Endobronchial Robertshaw Left	39 Fr	1
1-7370-41	Endobronchial Robertshaw Left	41 Fr	1

DUAL LUMEN (RIGHT)

ITEM	DESCRIPTION	SIZE	PK
1-7371-35	Endobronchial Robertshaw Right	35 Fr	1
1-7371-37	Endobronchial Robertshaw Right	37 Fr	1
1-7371-39	Endobronchial Robertshaw Right	39 Fr	1
1-7371-41	Endobronchial Robertshaw Right	41 Fr	1

LIT-10037 Rev.1