


INTRODUCER/BOUGIE

LATEX FREE | STERILE | SINGLE USE | DISPOSABLE | CE MARK


Latex free introducers by SunMed® are both firm and flexible for ease of insertion. They are calibrated for accurate entry depth and available in adult and pediatric sizes.



- Low density polyethylene provides proper stiffness for ease of insertion
- Calibrated to ensure accurate distance insertion
- Pediatric size features calibrated markings on both ends allowing it to be used as a straight or coudé tip
- Pediatric size fits in 4 mm to 6 mm tubes
- Adult size fits in 6 mm to 11 mm tubes

ENDOTRACHEAL TUBE INTRODUCER

ITEM	DESCRIPTION	SIZE	PK
9-0211-70	Pediatric, Coudé Tip	10 Fr x 70 cm	10
9-0212-70	Adult, Coudé Tip	15 Fr x 70 cm	10
9-0212-72	Adult, Straight Tip	15 Fr x 70 cm	10



MALLEABLE ENDOTRACHEAL TUBE INTRODUCER

- Combines benefits of an introducer and stylette into one unit
- Pediatric size features calibrated markings on both ends allowing it to be used as a straight or coudé tip

ITEM	DESCRIPTION	SIZE	PK
9-0213-82	Pediatric Malleable, Coudé & Straight Tip	10 Fr x 70 cm	10
9-0213-92	Adult, Malleable, Coudé Tip	15 Fr x 70 cm	10

BOUGIE-TO-GO™

LATEX FREE | STERILE | SINGLE USE | DISPOSABLE | CE MARK

The Bougie-To-Go™ latex-free introducer by SunMed® was developed with first responders in mind. Designed with the same firm and flexible quality as our traditional introducer, the Bougie-To-Go features compact packaging. It fits in paramedic bags, pockets and smaller storage spaces.



BOUGIE-TO-GO™ ENODOTRACHEAL INTRODUCER

- Innovative compact packaging designed specifically for use by EMS professionals
- Low density polyethylene provides proper stiffness for ease of insertion
- Calibrated to ensure accurate distance insertion
- Adult size fits in 6 mm to 11 mm tubes

ITEM	DESCRIPTION	SIZE	PK
9-0212-82	Adult, Coudé Tip	15 Fr x 60 cm	10



Unwraps without distorting introducer shape; keeps tip in correct insertion position