

# FIBER OPTIC GREENLINE® HANDLE

## LATEX FREE

Choose from a variety of handles for use with SunMed® brand fiber optic laryngoscope reusable blades. All are compatible with Green System blades and comply with ISO 7376 Green Standard.



## GREENLINE® HANDLE

- Chrome plated brass and surgical stainless steel handles are compatible with all "Green System" fiber optic blades.
- Bright Halogenated Xenon lamp

## CHROME PLATED


ITEM	DESCRIPTION	SIZE	BATTERIES	BOX
5-0236-09	Chrome Plated Handle	Medium	2 "C"	1
5-0236-10	Chrome Plated Handle	Penlite	2 "AA"	1
5-0236-11	Chrome Plated Handle	Stubby	2 "AA"	1
5-0236-15	Chrome Plated Handle	Large	2 "D"	1

## SURGICAL STAINLESS STEEL

ITEM	DESCRIPTION	SIZE	BATTERIES	BOX
5-0236-39	Stainless Steel Halogen Handle	Medium	2 "C"	1
5-0236-30	Stainless Steel Halogen Handle	Penlite	2 "AA"	1
5-0236-31	Stainless Steel Halogen Handle	Stubby	2 "AA"	1
5-0236-35	Stainless Steel Halogen Handle	Large	2 "D"	1
5-0236-38	Stainless Steel Halogen Handle	Mini	Lithium "123"	1
5-0236-36	Stainless Steel Halogen Handle	Micro	Lithium "123"	1

## PREMIUM LED HANDLE

- Premium LED handles have a lifetime warranty
- Low power use prolongs battery life
- Compliant with ISO standard 7376



ITEM	DESCRIPTION	SIZE	BATTERIES	BOX
5-0236-80	Premium LED Handle	Penlite	2 "AA"	1
5-0236-89	Premium LED Handle	Medium	2 "AA"	1

## SUNBRITE™ LED HANDLE

- Eliminates any heat in handle
- Cool, white illumination
- Low power use prolongs battery life



ITEM	DESCRIPTION	SIZE	BATTERIES	BOX
5-0236-90	SunBrite LED Handle	Penlite	2 "AA"	1
5-0236-91	SunBrite LED Handle	Stubby	2 "AA"	1
5-0236-99	SunBrite LED Handle	Medium	2 "C"	1

\*Batteries not included