

SunOne[™]
SYSTEM

GREENLINE®/D™ FIBER OPTIC LARYNGOSCOPES

LATEX FREE | STERILE | SINGLE PATIENT USE | DISPOSABLE

GREENLINE®/D™ GREEN HEEL BLADES

- Reduce the risk of cross-contamination with single patient use blades
- Constructed of non-painted 303/304 surgical stainless steel
- Acrylic stem facilitates maximum LED illumination for superior, clear visualization
- Soft matte finish prevents reflection and glare during tissue visualization
- Durable nylon blade heel with two ball bearing connections
- Extended safety heel design prevents blade-to-handle contact
- Fully disposable; product recycling available

5 YEAR STERILITY

Single patient use,
fully disposable fiber
optic laryngoscopes



 ASK ABOUT PRODUCT RECYCLING

GREENLINE®/D™ LED HANDLES

- Reduce the risk of cross-contamination with single patient use handles
- Textured surface ensures secure grip
- Mercury-free, pre-installed batteries
- Contained unit is fully disposable; product recycling available

5 YEAR STERILITY + SHELF LIFE



GREENLINE®/D™ GREEN HEEL BLADES



GREEN HEEL MACINTOSH

ITEM	DESCRIPTION	LENGTH	HEIGHT	SIZE
5-5332-01	Infant	92 mm	22 mm	1
5-5332-02	Child	100 mm	21 mm	2
5-5332-03	Medium Adult	130 mm	22 mm	3
5-5332-35	Ext. Med. Adult	144 mm	22 mm	3.5
5-5332-04	Large Adult	155 mm	25 mm	4

GREEN HEEL MILLER

ITEM	DESCRIPTION	LENGTH	HEIGHT	SIZE
5-5333-20	Premature	67 mm	10 mm	00
5-5333-00	Neonate	77 mm	10 mm	0
5-5333-01	Infant	102 mm	10 mm	1
5-5333-02	Child	153 mm	12 mm	2
5-5333-03	Medium Adult	195 mm	14 mm	3
5-5333-04	Large Adult	205 mm	18 mm	4

FLEX TIP BLADES

ITEM	DESCRIPTION	LENGTH	HEIGHT	SIZE
5-5382-03	Medium Adult	195 mm	18 mm	3
5-5382-04	Large Adult	205 mm	18 mm	4



ROBERTSHAW BLADES

ITEM	DESCRIPTION	LENGTH	HEIGHT	SIZE
5-5337-00	Neonate	94 mm	11 mm	0
5-5337-01	Infant	112 mm	13 mm	1

GREENLINE®/D™ LED HANDLES

ITEM	DESCRIPTION	STERILE	SIZE
5-5338-09	Textured Grip	Sterile	Medium
5-5338-10	Textured Grip	Sterile	Penlite
5-5338-11	Textured Grip	Sterile	Stubby
5-5338-12	Textured Grip	Sterile	Standard
5-0236-69	Ribbed Grip	Non-Sterile	Large
5-0236-70	Smooth Grip	Non-Sterile	Penlite



LIT-10051 Rev.1