


GREENLINE® FIBER OPTIC AMERICAN BLADE

LATEX FREE | REUSABLE

The SunMed™ Fiber Optic American profile laryngoscope blades have a soft matte finish that virtually eliminates reflection and glare. The fiber optic light bundle provides bright illumination and the blades are made of 303/304 surgical stainless steel, not chrome plated. They are Green System compatible and ISO standard 7376 compliant.



MACINTOSH FIBER OPTIC AMERICAN PROFILE

- The flange begins at the lock of the blade, continues to where the curvature changes direction, and tapers off
- Offers a wide variety of sizes to accommodate any patient need

ITEM	DESCRIPTION	SIZE	A	B	PK
5-5231-00	Newborn	0	80 mm	15 mm	1
5-5231-01	Infant	1	87 mm	15 mm	1
5-5231-02	Child	2	108 mm	15 mm	1
5-5231-03	Medium Adult	3	128 mm	24 mm	1
5-5231-35	Ext. Med. Adult	3.5	144 mm	22 mm	1
5-5231-04	Large Adult	4	159 mm	15 mm	1

MILLER FIBER OPTIC AMERICAN PROFILE

- The blade's distal end curve improves lifting of the epiglottis
- Designed to facilitate greater exposure of larynx
- Light-on-right profile

ITEM	DESCRIPTION	SIZE	A	B	PK
5-5230-02	Child	2	155 mm	13 mm	1
5-5230-03	Medium Adult	3	195 mm	13 mm	1

