

# GREENLINE® FIBER OPTIC ENGLISH BLADE

LATEX FREE | REUSABLE



The SunMed™ Fiber Optic English profile laryngoscope blades have a soft matte finish that virtually eliminates reflection and glare. The fiber optic light bundle provides bright illumination. They are made of 303/304 surgical stainless steel, not chrome plated, are Green System compatible and ISO standard 7376 compliant.



## MACINTOSH FIBER OPTIC ENGLISH PROFILE

- The flange begins at the lock of the blade and continues to the tip
- Offers a wide variety of sizes to accommodate any patient need, including a left handed option

ITEM	DESCRIPTION	SIZE	A	B	PK
5-5232-01	Infant	1	92 mm	22 mm	1
5-5232-02	Child	2	100 mm	21 mm	1
5-5232-03	Medium Adult	3	130 mm	22 mm	1
5-5232-04	Medium Adult	4	155 mm	25 mm	1
5-5232-05	Extra Lg. Adult	5	175 mm	22 mm	1
5-5243-02	Left Hand Adult	3	130 mm	22 mm	1



## MILLER FIBER OPTIC ENGLISH PROFILE

- The blade's distal end curve improves lifting of the epiglottis
- Designed to facilitate greater exposure of larynx
- Reduced flange helps minimize trauma

ITEM	DESCRIPTION	SIZE	A	B	PK
5-5233-20	Premature	00	67 mm	10 mm	1
5-5233-00	Newborn	0	77 mm	10 mm	1
5-5233-01	Infant	1	102 mm	10 mm	1
5-5233-15	Small Child	1.5	130 mm	12 mm	1
5-5233-02	Child	2	153 mm	12 mm	1
5-5233-03	Medium Adult	3	195 mm	12 mm	1
5-5233-04	Large Adult	4	205 mm	16 mm	1

