

# HEART® CONTINUOUS NEBULIZER

LATEX FREE | SINGLE PATIENT USE | DISPOSABLE | CE MARK

The family of HEART® Continuous Nebulizers offers 1 to 8 hours of treatment time. The continuous treatment can replace three separate small volume nebulizer treatments, improving workflow for the caregiver.



HEART Hi-Flo Nebulizer



IV HEART™ Nebulizer



MiniHEART-LoFlo® Nebulizer



MiniHEART-HiFlo® Nebulizer with Mask



UniHEART™ LoFlo Nebulizer

## HEART® CONTINUOUS NEBULIZERS

- The HEART high-output nebulizer provides up to 8 hours of therapy, flow rates of 10-15 LPM, and output up to 50 mL/hr
- IV HEART™ and UniHEART™ have an IV port for easy refills.
- UniHEART™ and MiniHEART-HiFlo® are ideal for emergency department.
- MiniHEART-LoFlo® designed for ventilated patients, and has an IV port for easy refills.

ITEM	THERAPY TIME	RESERVOIR	PARTICLE SIZE
HEART	8 hours (10 LPM)	240 mL	2-3 µm
IV HEART	4 hours (10 LPM)	100 mL	2-3 µm
MiniHEART-LoFlo	3 hours (2 LPM)	30 mL	2-3 µm
MiniHEART-HiFlo	1 hour (8 LPM)	30 mL	2-3 µm
UniHEART	2 hours (2 LPM)	10 mL	2-3 µm

## HEART® Nebulizer

ITEM	DESCRIPTION	PK
100609	HEART	18
100610	HEART Direct Connect	18
100619	Adult HEART Direct Connect Set-Up Kit – HEART, Adult Mask, 8' 22mm Hose	18
100629	Pediatric HEART Direct Connect Set-Up Kit – HEART, Pediatric Mask, 8' 22mm Hose	18
100639	Adult HEART Set-Up Kit – HEART, Adult Mask, 8' 22mm Hose & 7' Oxygen Tubing with Nut Connector	18
100649	Pediatric HEART Set-Up Kit – HEART, Pediatric Mask, 8' 22mm Hose & 7' Oxygen Tubing with Nut Connector	18

## IV HEART™ Nebulizer

ITEM	DESCRIPTION	PK
101020	IV HEART	25
101021	Adult IV HEART Set-Up Kit – IV HEART, Adult Mask, 8' 22mm Hose & 7' Oxygen Tubing with Nut Connector	18

## MiniHEART-LoFlo® Nebulizer

ITEM	DESCRIPTION	PK
100611	MiniHEART-LoFlo with IV Drip Port	25
100611A	MiniHEART-LoFlo with IV Drip Port, with Tee 6" Flex Hose & 7' Oxygen Tubing with Nut Connector	25
100611B	MiniHEART-LoFlo	25
100611C	MiniHEART-LoFlo with extension Set & Neonatal Valved Tee	25
100623	MiniHEART-LoFlo with extension Set & Adult Valved Tee	25

## MiniHEART-HiFlo® Nebulizer

ITEM	DESCRIPTION	PK
100612	MiniHEART-HiFlo Continuous Nebulizer	25
100612A	MiniHEART-HiFlo Continuous Nebulizer with Adult Mask & 7' Oxygen Tubing with Nut Connector	25
100613A	MiniHEART-HiFlo Continuous Nebulizer with Pediatric Mask & 7'Tubing	25
100614A	MiniHEART-HiFlo Continuous Nebulizer with Super Spike Mask & 7'Tubing	25

## UniHEART™ Low-Flow Nebulizer

ITEM	DESCRIPTION	PK
100850	UniHEART Lo-Flow	50
100940	UniHEART Lo-Flow with Adult Mask & 7' Oxygen Tubing with Nut Connector	25
100942	UniHEART Lo-Flow with Pediatric Super Spike Mask & 7'Tubing	25