


WESTMED™ ETCO₂ CANNULAS WITH REFLECTIVE CONNECTORS

LATEX FREE | SINGLE PATIENT USE | DISPOSABLE | DEHP FREE

SunMed® offers the Westmed™ EtCO₂® cannula line, a comprehensive range of adult and pediatric capnography solutions. The nasal interface is bifurcated to effectively deliver O₂ while capturing CO₂ levels. The specially designed oral/nasal cannula interface supports patients who breathe through their mouths.



WESTMED™ ETCO₂ CANNULAS WITH REFLECTIVE CONNECTORS

- Advanced design helps minimize dead space in the sampling line, resulting in a more accurate capnogram waveform.
- Comfort Soft Plus® tubing helps reduce Hospital-Acquired Pressure Ulcers (HAPU), while promoting patient wellness. (Clinical studies available)^{1, 2}
- The unique filtration system helps to prevent moisture from entering the line and infiltrating the monitor.
- Multiple models, designs and custom configurations available:
 - Bifurcated O₂ and CO₂ prongs
 - Oral/Nasal design for mouth-breathers
 - 7', 10', and 14' sampling line lengths



Compatible with Microstream™, CapnoStream™, and Oridion™ capnography monitors

1. Mullen J, Johnston J, Bryner L, Pearson C, Weideman Y. A multidisciplinary approach to the prevention of ear pressure ulcers. JWON 41. S54-S55
 2. Lukose M, Thomas L, Camacho-Del Rio G, et al. No more sore ears: Nasal cannula related pressure ulcers down to zero. WOCN 49th Annual Conference. 2017.

0913 Oral/Nasal Cannula

WESTMED™ ETCO₂ CANNULAS WITH REFLECTIVE CONNECTORS

ITEM	DESCRIPTION	SIZE	CO ₂ LENGTH	O ₂ LENGTH	PK
0912	Nasal Cannula	Adult	10'	10'	25
0913	Oral/Nasal Cannula	Adult	10'	10'	25
0914	Nasal Cannula	Adult	14'	14'	25
0915	Oral/Nasal Cannula	Adult	14'	14'	25
0916	Nasal Cannula	Pediatric	10'	10'	25
0917	Oral/Nasal Cannula	Adult	7'	7'	25
0918	Oral/Nasal Cannula	Adult	2"	7'	25
0920	Nasal Cannula with Gas Sampling Tee	Adult	7'	7'	25
0921	Oral/Nasal Cannula with Gas Sampling Tee	Adult	7'	7'	25
0930	CO ₂ /O ₂ Mask with 90° Swivel Connector	Adult	10'	10'	25
0931	CO ₂ /O ₂ Mask with 90° Swivel Connector	Pediatric	10'	10'	25

GAS SAMPLING LINE WITH 90° SWIVEL CONNECTOR

ITEM	DESCRIPTION	LENGTH	PK
0940	GSL	10'	25
0942	GSL, Gas Sampling Tee	10'	25