


47SOFT DIVIDED ET_{CO}₂ SAMPLING CANNULA

LATEX FREE | SINGLE PATIENT USE | DISPOSABLE | DEHP FREE

The Salter Labs® divided capnography cannulas are now available with soft tubing. Salter Labs design provides true EtCO₂ sampling and proven O₂ delivery without limitations at higher O₂ rates. The divided design separates gases and Salter Eyes™ maintain sampling even during a cannula occlusion. Available in a multitude of configurations for Luer style connectors.



47SOFT DIVIDED ET_{CO}₂ SAMPLING CANNULA

- Salter Labs® soft headset tubing has been shown to improve patient comfort². Pliable, stretchable tubing curves around the cheeks and ears.
- 3-Channel tube design to resist kinking and reduce the risk of oxygen or EtCO₂ flow loss.
- SalterStyle™ Divided Nasal Sampling Cannulas provide true EtCO₂ sampling and effective O₂ delivery without limitations¹
- Salter Eyes™ provide safety apertures in each prong to reduce the possibility of occlusions.³

References: 1. Ebert TJ, Novallja J, Uhrich TD, Barney JA. The effectiveness of oxygen delivery and reliability of carbon dioxide waveforms: a crossover comparison of 4 nasal cannulae. *Anesth Analg*. 2015;120(2):342-438. 2. Salter Labs Data on File. 3. Salter Labs Technical Report 006.

47SOFT DIVIDED ETCO₂ SAMPLING CANNULA

SOFT DIVIDED ETCO₂ SAMPLING CANNULA

ITEM	CO ₂ CONNECTION	O ₂ CONNECTION	CO ₂ LINE	O ₂ LINE	SIZE	PK
47SOFT-7-7-25	Male Luer Style CO ₂	Standard	7'	7'	Adult	25
47SOFT-13-13-25	Male Luer Style CO ₂	Standard	13'	13'	Adult	25
47SOFT-7-0-25	Male Luer Style CO ₂	Standard	2"	7'	Adult	25
47SOFT-TG-7-7-25	Male Luer Style CO ₂	Thread Grip	7'	7'	Adult	25
47SOFT-TG-14-14-25	Male Luer Style CO ₂	Thread Grip	14'	14'	Adult	25
47SOFT-TG-7-0-25	Male Luer Style CO ₂	Thread Grip	2"	7'	Adult	25
47SOFT-PED-7-0-25	Male Luer Style CO ₂	Standard	2"	7'	Pediatric	25
47SOFT-PED-7-2-25	Male Luer Style CO ₂	Standard	2"	7'	Pediatric	25
47SOFT-PED-7-7-25	Male Luer Style CO ₂	Standard	7'	7'	Pediatric	25
47SOFT-TG-PED-7-7-25	Male Luer Style CO ₂	Thread Grip	7'	7'	Pediatric	25
47SOFT-F-7-7-25	Female Luer Style CO ₂	Standard	7'	7'	Adult	25
47SOFT-F-14-14-25	Female Luer Style CO ₂	Standard	14'	14'	Adult	25
47SOFT-F-7-0-25	Female Luer Style CO ₂	Standard	2"	7'	Adult	25
47SOFT-FTG-7-7-25	Female Luer Style CO ₂	Thread Grip	7'	7'	Adult	25
47SOFT-FTG-7-0-25	Female Luer Style CO ₂	Thread Grip	2"	7'	Adult	25
47SOFT-FTG-PED-7-7-25	Female Luer Style CO ₂	Thread Grip	7'	7'	Pediatric	25