

ARTERIAL BLOOD GAS (ABG) SAMPLING

LATEX FREE | SINGLE USE | DISPOSABLE

Arterial Blood Gas (ABG) sampling is an invasive procedure that can result in accidental needlestick injury and exposure to bloodborne pathogens. The Pulset™ ABG Syringe features SafetyTip™ Needlestick Protection to help prevent accidental needlestick injuries in the care setting. The SafetyTip is engineered as an integral part of the needle to comply with all FDA, OSHA and NIOSH guidelines. The Pulset has other additional features making it one of the easiest and quickest ABG syringes to use in healthcare today.



Puncture Kit



A-Line Kit



PULSET™ ADVANTAGES

SHARPS PROTECTION DEVICES (SAFETYTIP™)

- The SafetyTip comes attached to the needle with bevel in the ready/up position
- The SafetyTip device folds back along the syringe for optimal visualization of the insertion site
- Single hand activation

VENTED PLUNGER

- Eliminates the need for manual clearance of air bubbles
- Allows for fewer steps between needle insertion and sending the sample to the lab

DRY BALANCED LITHIUM-ZINC HEPARIN

- Eliminates calcium electrolyte errors
- Easily mixes with blood to prevent clotting

THIN-WALLED NEEDLE TECHNOLOGY

- Flows as quickly as standard needles
- Smaller puncture site means less pain for the patient

Distinctive by design, SafetyTip™ provides optimum visualization of the insertion site. With the needle bevel facing up, the protective sheath aligns parallel to the syringe barrel, resulting in a clear, unobstructed view. In addition, the patented Dual Locking design promotes your safety by securely covering the needle with an audible click.

1CC PULSET SYRINGE, VENTED

ITEM	DESCRIPTION	ALCOHOL PREP	IODINE PREP	GAUZE	ICE BAG	BANDAGE	Q-CORK	PATIENT LABEL	QTY
3103-91*	23g x 1" Needle with SafetyTip, 25u BH								100
3105-91*	25g x 5/8" Needle with SafetyTip, 25u BH								100
3165-91	Modified Kit, 25g x 5/8" Needle with SafetyTip, 25u BH	X	X	X	X	X		X	100
3173-91	Full Kit: 25g x 5/8" Needle 23g x 1" Sterile Unattached Needle, 25u BH	X	X	X	X	X	X	X	100

3CC PULSET SYRINGE, VENTED

ITEM	DESCRIPTION	ALCOHOL PREP	IODINE PREP	GAUZE	ICE BAG	BANDAGE	Q-CORK	PATIENT LABEL	QTY
3302-91*	22g x 1" Needle with Safety Tip, 25u BH								100
3303-91*	23g x 1" Needle with Safety Tip, 25u BH								100
3362-91	Modified Kit, 22g x 1" Needle with SafetyTip, 25u BH	X	X	X	X	X		X	100
3363-91	Modified Kit, 23g x 1" Needle with SafetyTip, 25u BH	X	X	X	X	X		X	100
3373-91	Full Kit: 23g x 1" Needle with SafetyTip, 22g x 1.5" Sterile Unattached Needle with Safety Tip, 25u BH	X	X	X	X	X		X	100

CORD GAS KIT

ITEM	DESCRIPTION	ALCOHOL PREP	IODINE PREP	GAUZE	ICE BAG	BANDAGE	Q-CORK	PATIENT LABEL	QTY
3182-91	Cord Gas Kit, (2) 22g x 1" Syringes with Needles with Safety Tip, 25u BH				X			X (2)	100

A-LINE WITH BLUNT TIP CANNULA

ITEM	DESCRIPTION	QTY
3100-30*	A-Line 1cc Syringe Blunt Tip Cannula 25u BH	100
3100-31*	A-Line/Pulset 1cc ABG, Blunt Tip Cannula 25u BH	100
3300-30*	A-Line 3cc Syringe Blunt Tip Cannula 25u BH	100
3300-31*	A-Line/Pulset 3cc ABG, Blunt Tip Cannula 25u BH	100

A-LINE WITH NO NEEDLE/CANNULA

ITEM	DESCRIPTION	QTY
3100-20*	A-Line 1cc Syringe 3cc Waste Syringe 25u BH	100
3100-25*	A-Line 1cc Syringe 25u BH	100
3300-21*	A-Line 3cc Syringe 5cc-6cc Waste Syringe 25u BH	100
3300-25*	A-Line 3cc Syringe 25u BH	100

*Items include CE Mark.