


EtCO₂ SalterSTAT™ Overview

EtCO₂ monitoring is a noninvasive breath by breath analysis of a patient's ventilatory status. Capnography is rapidly expanding into intensive-care and procedural care settings (GI Suites) to provide safer patient monitoring.

As the leader in EtCO₂ consumables utilizing standard luer connectors, we now offer cost-effective solutions for monitors requiring reflective connectors.

SalterSTAT™ Reflective Series

For Monitors Requiring Reflective Connectors



Compatible Reflective Connector

- **Reputation:** Created by a trusted source for respiratory care products for over 40 years.
- **Compatibility:** Full range of configurations designed to be compatible with monitors that require reflective connectors including: Philips™, Covidien™, Medtronic™, Physio-Control™, Mindray™, Zoll™, Datascope™, HP, Respironics™
- **Flexibility:** No annual contract required
- **Versatility:** Compatible sampling cannulas, lines and masks
- **Ergonomic Design:** Provides superior comfort and fit on both adult and pediatric faces
- **Moisture Management:** Standard in-line 0.2 micron hydrophobic filter to protect sensitive equipment. Also available with Nafion tubing for long term monitoring

TOP USAGE PRODUCT CODES – ACUTE CARE AND EMS Competitive Cross Matrix

ITEM	DESCRIPTION	MEDTRONIC	PHILIPS	PHYSIO-CONTROL	ZOLL
4MSF1-7-6-25	Oral-Nasal sampling with 7' O ₂ , 6' CO ₂ line	009822	989803129751	11996-000163	8300-000528-01
4MSF3-7-6-25	Nasal sampling with 7' O ₂ , 6' CO ₂ line	006912	989803179101	–	–
4MSF3-6-25	Nasal sampling 6' CO ₂ line	008174	–	–	–
4MSF1-7-13-25	Oral-Nasal sampling with 7' O ₂ , 13' CO ₂ line	009826	989803129781	11996-000165	8300-0524-01
4MSF5-13-25	Sampling line with ETT adapter, 13' CO ₂ line	007768	–	–	8300-0520-01

TOP USAGE PRODUCT CODES – CHILDREN'S HOSPITAL Competitive Cross Matrix

ITEM	DESCRIPTION	MEDTRONIC	PHILIPS	PHYSIO-CONTROL	ZOLL
4MSF2-PED-7-6-25	Oral-Nasal sampling with 7' O ₂ , 6' CO ₂ line	007269	989803129761 M2524A	–	–
4MSF3-6-25	Nasal sampling 6' CO ₂ line	008175	–	–	–
4MSF2-L-PED-7-6-25	Oral-Nasal sampling with 7' O ₂ , 6' CO ₂ line	010582	989803177971	–	–
4MSF5-13-25	Sampling line with ETT adapter, 13' CO ₂ line	007768	–	–	8300-0520-01


More configurations available, view our STATOrange Product selection guide.

For Long Term Use

Long term use EtCO₂ cannula options include Nafion Tubing for moisture management.



Long Term Use
Oral-Nasal Cannula



Long Term Use
Nasal Cannula

LONG TERM


ITEM	DESCRIPTION	MEDTRONIC	PHILIPS	PHYSIO-CONTROL	ZOLL
4MSF1-L-7-6-25	Long term use Oral-Nasal sampling with 7' O ₂ , 6' CO ₂ line	010433	989803177951	-	-
4MSF2-L-PED-7-6-25	Long term use Oral-Nasal sampling with 7' O ₂ , 6' CO ₂ line,	010582	989803177971	-	-
4MSF5-L-6-25	Long term use sampling line with ETT adapter, 6' CO ₂ line	010787	-	-	-
4MSF8-L-13-25	Long term use sampling line with male luer-lock, 13' CO ₂ line	-	-	-	-

More configurations available, view our STATOrange Product selection guide.

EtCO₂ SalterSTAT™ Overview

With the recent Salter Labs® acquisition, SunMed® is now in a position to offer an extensive EtCO₂ cannula technology with industry demanded performance and competitive pricing as our feedback from clinicians state.

Performance Comparison



Waveform SalterSTAT Reflective Series

Waveform Microstream

THE SALTERSTAT™ REFLECTIVE SERIES RECEIVED EXCELLENT CLINICAL EVALUATION FEEDBACK SUCH AS:

- Provides clean, crisp waveforms
- Reliable, accurate values throughout procedure
- Exceptional comfort
- Device performs well with NRB mask
- Performs well during CPR

