

AirFlow™



 SunMed®

Pictured Adult AF5140MB

The SunMed® AirFlow™ Resuscitation Bag offers new, enhanced features for higher levels of accuracy and safety.

- Redesigned manometer to improve visibility and enhance accuracy.
- Improved accuracy on our patented, integrated color-coded pressure manometer.¹
- MRI Conditional up to 7 tesla.
- Smaller patient valve allows for decreased dead space < 8 mL.²
- Sampling port allows for EtCO₂ monitoring.

1. ± 3 cm H₂O Up to 15cm H₂O
2. < 5mL + 10% of delivered volume